


Product Datasheet

6x-His Tag antibody (orb108772)

Description	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Detects only 6x-His tag by Western blot. Specific for His tags no longer than 6 histidines. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG2a
Clonality	Monoclonal
Source	Hybridoma produced by the fusion of splenocytes from mice immunized with His-His-His-His-His-His synthetic peptide and mouse myeloma cells.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of His tagged prot...