

Product Datasheet

Mouse anti-Interleukin 6, clone 10C12 (monoclonal) (orb1087729)

Description	Mouse anti-Interleukin 6, clone 10C12 (monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tissue culture supernatant with 15mM Sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	IHC-P (1:50)
Expiration Date	12 months from date of receipt.