

Product Datasheet

Mouse anti-Helicobacter pylori, Clone ULC3R (monoclonal) (orb1087733)

Description	Mouse anti-Helicobacter pylori, Clone ULC3R (monoclonal)
Species/Host	Mouse
Reactivity	Bacteria, Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tissue culture supernatant with 15mM Sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.