

Product Datasheet

Lysophosphatidic Acid Receptor 3 CT antibody (orb108774)

Description

Rabbit polyclonal to Lysophosphatidic Acid

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Antibody raised against a unique synthetic peptide derived from the C-terminal of the human EDG7 protein.

Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Detects recombinant EDG-7 receptors by Western blot. Due to low expression of EDG receptors, we recommend use of Pierce Femto Signal substrate for western blot development. Optimal concentration should be evaluated by serial dilutions.

Isotype

IgG

Uniprot ID
Q9UBY5
Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.


EDG-7
(40 kDa)

Western
blot
analysis of
RH7777
cell lys...