

Product Datasheet

Mouse anti-Cytokeratin 5 (Liquid Primary Antibody), clone XM26 (monoclonal) (orb1087760)

Description	Mouse anti-Cytokeratin 5 (Liquid Primary Antibody), clone XM26 (monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tissue culture supernatant with Sodium azide
Concentration	Greater than or equal to 21 mg/L
Storage	2-8°C
Note	For research use only
Isotype	IgG1, kappa
Clonality	Monoclonal
Dilution Range	IHC-P (1:100)
Expiration Date	12 months from date of receipt.