

Product Datasheet

Rit antibody (orb108777)

Description	Rabbit polyclonal to Rit
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	His-tagged recombinant human Rit protein
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody detects recombinant Rit protein expressed in HEK-293 cells by Western blot analysis. Shows minimal cross-reactivity with Rin protein and other Ras proteins. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Uniprot ID	Q92963
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.