

Product Datasheet

Mouse anti-Beta-Dystroglycan, clone 43DAG1/8D5 (monoclonal) (orb1087799)

Description	Mouse anti-Beta-Dystroglycan, clone 43DAG1/8D5 (monoclonal)
Species/Host	Mouse
Reactivity	Amphibian, Canine, Gallus, Hamster, Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Lyophilized
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.

