



Product Datasheet

Interferon gamma antibody (orb10878)


Description	Interferon gamma antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Interferon gamma. Please contact us for the exact immunogen sequence. The peptide is available as orb393916.
Target	Interferon gamma
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	19 kDa
Uniprot ID	P01581, P01580
Entrez	15978
Dilution Range	IF/ICC: 1:100-600, WB: 1:200-1000, IHC-P: 1:100-600 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat test...



IHC-P staining of mouse kidney tissue
US...