

## Product Datasheet

### Proximal Nephrogenic Renal Antigen, cIn PN-15 antibody (orb108788)

**Description**

Mouse monoclonal to Proximal

15

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

IHC-P

**Applications**
**Preservatives**

Provided as solution in phosphate buffered saline with 0.08% sodium azide

**Form/Appearance**

Provided as solution in phosphate buffered saline with 0.08% sodium azide

**Concentration**

See vial for concentration

**Storage**

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

**Note**

For research use only

**Application notes**

This antibody can be used for the differential diagnosis of carcinomas by Western blotting and immunohistochemistry on frozen and paraffin-embedded tissues. Protease digestion prior to antibody incubation enhances staining of tissue sections. Optimal concentration should be evaluated by serial dilutions.

**Isotype**

IgG2b

**Clonality**

Monoclonal

**Source**

Hybridoma produced by the fusion of splenocytes from BALB/c mice immunized with renal cortical tissue extract and mouse myeloma cells

**Uniprot ID**

**Q16790**

**Hazard**
**Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.