

Product Datasheet

Keratin 13 & 10 antibody (orb108791)

Epidermal Reagents.

Mouse monoclonal to Keratin 13 & 10

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Hybridoma produced by the fusion of splenocytes from mice immunized with cytoskeleton preparation from HeLa cells and mouse myeloma cells.
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody detects cytokeration protein numbers 13 and 10 (molecular weights 53 and 56.6 kDa, respectively, according to Moll's classification) by Western blot. Does not react with other cytokeratins or other proteins. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P13646
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.