

Product Datasheet

alpha-Keratin 5 & 8 antibody (orb108795)

Description

Mouse monoclonal to alpha-Keratin 5 & 8

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-Fr, WB

Applications

Immunogen

Hybridoma produced by the fusion of splenocytes from mice immunized with cytoskeleton preparation of HeLa cells and mouse myeloma cells.

Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Detects Keratin 5 & 8 by Western blot. Detects a 58 kDa band, corresponding to Keratin 5, and a 52.5 kDa band, corresponding to Keratin 8. Optimal concentration should be evaluated by serial dilutions.

Isotype

IgG1

Clonality

Monoclonal

Uniprot ID

P05787

Hazard

Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.