

Product Datasheet

Mouse anti-CD45RO, clone UCHL-1 (orb1087993)

Description	Mouse anti-CD45RO, clone UCHL-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2a-k
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.