

Product Datasheet

TCC, Human, mAb, clone aE11 (Monoclonal) Antibody (orb1088073)

Description	TCC, Human, mAb, clone aE11 (Monoclonal)
Species/Host	Mouse
Reactivity	Equine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FC, IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.