

Product Datasheet

Mouse anti-Cytokeratin 7, Keratin K7 (Monoclonal) (orb1088093)

Description	Mouse anti-Cytokeratin 7, Keratin K7 (Monoclonal)
Species/Host	Mouse
Reactivity	Canine, Feline, Goat, Hamster, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,,,
Expiration Date	12 months from date of receipt.