

Product Datasheet

Mouse anti-PD-L1, clone 405.9A11 (Monoclonal) (orb1088111)

Description	Mouse anti-PD-L1, clone 405.9A11 (Monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tissue culture supernatant with sodium azide
Concentration	lot specific
Storage	2-8°C
Note	For research use only
Isotype	IgG1k
Clonality	Monoclonal
Dilution Range	IHC-P (1:50-1:100)
Expiration Date	12 months from date of receipt.