

Product Datasheet

Human DPP3(Dipeptidyl peptidase 3) ELISA Kit (orb1088145)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, saliva, cell lysates, follicular fluid and other biological fluids
Uniprot ID	Q9NY33
Sensitivity	0.061 ng/mL
Expiration Date	Please enquire.

