

Product Datasheet

Rat sFLT(Soluble Fms-Like Tyrosine Kinase Receptor 1) ELISA Kit (orb1088159)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	78.13-5000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, saliva, cell lysates, follicular fluid and other biological fluids
Sensitivity	29 pg/mL
Expiration Date	Please enquire.

