

Product Datasheet

Mouse CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit (orb1088178)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P39689
Sensitivity	0.059 ng/mL
Expiration Date	Please enquire.

