

Product Datasheet

Zebrafish BDNF(Brain Derived Neurotrophic Factor) ELISA Kit (orb1088185)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Zebrafish
Range	7.82-500 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Sensitivity	3.1 pg/mL
Expiration Date	Please enquire.

