

Product Datasheet

Dog TNFRSF5(Tumor Necrosis Factor Receptor Superfamily, Member 5) ELISA Kit (orb1088197)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Canine
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	Serum, plasma and other biological fluids.
Uniprot ID	Q7YRL5
Sensitivity	5.4 pg/mL
Expiration Date	Please enquire.

