



Product Datasheet

IGF1 Receptor antibody (orb10882)

Description	IGF1 Receptor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human IGF-1 receptor beta chain (1102-1204/1204 aa)
Target	IGF1R
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150 kDa
Uniprot ID	P08069
Dilution Range	Flow-Cyt=1µg/Test, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P of mouse spleen tissue (dilution a...