

Product Datasheet

TGFbeta receptor III antibody (orb108833)

Description

Sheep polyclonal to TGFbeta receptor III

Species/Host

Sheep

Reactivity

Gallus

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

Recombinant protein derived from the extracellular domain of the chicken TGF III receptor protein. Protein shares significant homology with the mouse, rat and human sequences.

Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Antibody can be used for Western blotting (5-10 µg/ml), immunoprecipitation (5 µg/mg of protein lysate and immunohistochemistry on paraffin embedded tissue sections (5-10 µg/ml). Optimal concentration should be evaluated by serial dilutions.

Isotype

IgG

Uniprot ID


P16047

Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of extraemb...