

Product Datasheet

Prostaglandin-E2 receptor EP3 antibody (orb108837)

Description

Mouse monoclonal to Prostaglandin-E2 receptor

Species/Host

Mouse

Reactivity

Bovine, Human, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Western blot analysis of bovine brain ly...

Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

This antibody can be used for Western blot analysis (1-5 µg/ml). Optimal concentration should be evaluated by serial dilutions.

Isotype

IgG2a

Clonality

Monoclonal

Source

Hybridoma produced by the fusion of splenocytes from mice immunized with recombinant human EP3 receptor protein and mouse myeloma cells.

Uniprot ID

P34980

Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.



← EP3
(55 kDa)