



Product Datasheet

Anti-CD235a Biotin (orb1088488)

Description	CD235a (glycophorin A, GPA) is a transmembrane sialoglycoprotein expressed on erythrocytes and their...
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	FC, IHC-Fr, IHC-P
Immunogen	Membrane preparation from splenic hairy cell leukemia
Target	CD235a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P02724
Entrez	2993
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD235a positive eryt...



Flow cytometry surface staining pattern ...