



Product Datasheet

CD173 Antibody (orb1088490)

Description	CD173 (blood group antigen H2) is a fucosylated saccharide (Fuc-alpha-1-2-Gal-beta-1-4-GlcNAc-beta) ...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human thrombocytes
Target	CD173
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Dilution Range	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD173 positive eryth...



Flow cytometry surface staining pattern ...