



Product Datasheet

CD307a Antibody (orb1088507)


Description	CD307a is a type I transmembrane glycoprotein of the Fc receptor family. It contains two ITAM motifs...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	DNA-immunization followed by a boost with the CD307a transfected cells
Target	CD307a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q96LA6
Entrez	115350
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD307a
positive
lymp...



Flow
cytometry
surface
staining
pattern ...