

## Product Datasheet

### Recombinant Human ADP Protein (orb1088544)

|                            |                                                                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Recombinant Human adiponectin (ADP) protein is produced by Mammalian expression system and the targe...                                     |
| <b>Tested Applications</b> | CLIA, EIA, LF                                                                                                                               |
| <b>Target</b>              | ADP                                                                                                                                         |
| <b>Preservatives</b>       | PBS, pH 7.4.                                                                                                                                |
| <b>Form/Appearance</b>     | Supplied as a 0.2 µm filtered solution of PBS, pH 7.4                                                                                       |
| <b>Concentration</b>       | Lot specific(±20%)                                                                                                                          |
| <b>Storage</b>             | Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.                                     |
| <b>Tag</b>                 | His-tag                                                                                                                                     |
| <b>Note</b>                | For research use only                                                                                                                       |
| <b>Application notes</b>   | Use: Calibration/quality control. Quality control: Mol Mass: 26.4 kDa AP Mol Mass: 28.5 kDa Purity: ≥90% as determined by reducing SDS-PAGE |
| <b>Purity</b>              | > 95% as determined by SDS-PAGE.                                                                                                            |
| <b>Source</b>              | HEK293                                                                                                                                      |
| <b>MW</b>                  | 30 kDa                                                                                                                                      |
| <b>Uniprot ID</b>          | <b>Q15848</b>                                                                                                                               |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                              |