

Product Datasheet

Recombinant Human PCT Protein (orb1088570)

Description	Recombinant Human procalcitonin (PCT) protein is produced by E.coli expression system and the target...
Tested Applications	CLIA, EIA, LF
Target	PCT
Preservatives	PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4
Concentration	Lot specific(±20%)
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Tag	His-tag
Note	For research use only
Application notes	Use: Calibration/quality control. Quality control: Mol Mass: 14.0 kDa AP Mol Mass: 15.0 kDa Purity: ≥90% as determined by reducing SDS-PAGE
Purity	> 90% as determined by SDS-PAGE.
Source	E.coli
MW	12.8 kDa
Expiration Date	6 months from date of receipt.