

Product Datasheet

Podocin polyclonal antibody (KU) (orb1088651)

Description	Podocin polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding to Human Podocin.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 41 kDa
Uniprot ID	Q9NP85
Dilution Range	WB: 1:500~1:1000
Expiration Date	12 months from date of receipt.