

Product Datasheet

CHMP4B+CHMP4C Polyclonal Antibody (orb1088662)

Description	CHMP4B+CHMP4C Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CHMP4B/CHMP4C:101-200/224
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	25kDa
Uniprot ID	Q96CF2, Q9H444
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.