

Product Datasheet

Sphingosine 1 Phosphate 4 Receptor antibody (orb108867)

Description

Mouse monoclonal to Sphingosine 1 Phosphate 4

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, WB

Applications
Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

This antibody can be used for Western blotting (5-10 µg/ml), Immunofluorescence and Flow Cytometry. Optimal concentration should be evaluated by serial dilutions.

Isotype

IgG1

Clonality

Monoclonal

Source

Hybridoma produced by the fusion of splenocytes from mice immunized with Unique synthetic peptide derived from the C-terminal of the EDG6 (S1P4) protein and mouse myeloma cells.

Uniprot ID

095977 (Human)

Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.


Expiration Date

12 months from date of receipt.

98

31
19
17
11
6

Western blot analysis of cell lysate using...



IF analysis of CHO cells using Sphingosi...