

Product Datasheet

PLXNC1/CD232 Polyclonal Antibody (orb1088675)

Description	PLXNC1/CD232 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human PLXNC1/CD232:681-780/1568
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	172kDa
Uniprot ID	O60486
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.