

Product Datasheet

Wilms Tumor Protein Polyclonal Antibody (orb1088681)

Description	Wilms Tumor Protein Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Wilms Tumor Protein:301-400/449
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	55kDa
Uniprot ID	P19544
Dilution Range	IHC-P=1:100-500 Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.