

Product Datasheet

phospho-MCK10 (Tyr513) Polyclonal Antibody (orb1088706)

Description	phospho-MCK10 (Tyr513) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Equine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MCK10 around the phosphorylation site of Tyr513:PA(p-Y)RL
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	120 kD
Uniprot ID	Q08345
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.