

Product Datasheet

S100A7 STARD12 Polyclonal Antibody (orb1088708)

Description	S100A7 STARD12 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Equine, Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human S100A7:41-101/101
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	11kDa
Uniprot ID	P31151
Dilution Range	Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.