

Product Datasheet

RAB9 Polyclonal Antibody (orb1088711)

Description	RAB9 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human RAB9:65-170/201
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	23 kD
Uniprot ID	P51151
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.