

Product Datasheet

CD109 Polyclonal Antibody (orb1088714)

Description	CD109 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD109:901-1000/1445
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	154kDa
Uniprot ID	Q6YHK3
Dilution Range	Flow-Cyt=3ug/Test
Expiration Date	12 months from date of receipt.