

Product Datasheet

phospho-CaMK2 alpha (Thr305) Polyclonal Antibody (orb1088741)

Description	phospho-CaMK2 alpha (Thr305) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Canine, Gallus, Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CaMK2 alpha around the phosphorylation site of Thr305:IL(p-T)TM
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	54 kD
Uniprot ID	Q13554
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.