

Product Datasheet

SARS-CoV Nucleocapsid Polyclonal Antibody (orb1088752)

Description	SARS-CoV Nucleocapsid Polyclonal Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS N (nucleocapsid):1-42/422
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	46 kD
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.