

Product Datasheet

SARS-CoV-2 N protein Polyclonal Antibody (orb1088767)

Description	SARS-CoV-2 N protein Polyclonal Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -58°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~46 kDa
Uniprot ID	PODTC9
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.