

Product Datasheet

Nestin Polyclonal Antibody (orb1088794)

Description	Nestin Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse Nestin:1231-1330/1864
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -31°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~209kDa
Uniprot ID	Q6P5H2
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.