

Product Datasheet

MSGN1 Polyclonal Antibody (orb1088847)

Description	MSGN1 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Equine, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MSGN1:101-193/193
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -23°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~21 kDa
Uniprot ID	A6NI15
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.