

Product Datasheet

EAAT3 Polyclonal Antibody (orb1088850)

Description	EAAT3 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human EAAT3:451-524/524
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -26°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~58 kDa
Uniprot ID	P43005
Dilution Range	Flow-Cyt=1ug/Test not yet tested in other applications.
Expiration Date	12 months from date of receipt.