

Product Datasheet

IKB alpha Polyclonal Antibody (orb1088915)

Description	IKB alpha Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human NFKBIA:1-120/314
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	37 kD
Uniprot ID	Q9Z1E3
Dilution Range	WB=1:500-2000 IHC-P=1:100-500 Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.