

Product Datasheet

Ctip1 Polyclonal Antibody (orb1088961)

Description	Ctip1 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Gallus, Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Ctip1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	91 kD
Uniprot ID	Q9H165
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.