

Product Datasheet

Neurocan Polyclonal Antibody (orb1088997)

Description	Neurocan Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Neurocan
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	140kD
Uniprot ID	P55066
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.