

Product Datasheet

Phospho-PIM1 (Tyr309) Polyclonal Antibody (orb1089024)

Description	Phospho-PIM1 (Tyr309) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PIM1 around the phosphorylation site of Tyr309
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	45kD
Uniprot ID	P11309
Dilution Range	Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.