

Product Datasheet

Phospho-PERK (Thr980) Polyclonal Antibody (orb1089026)

Description	Phospho-PERK (Thr980) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from mouse PERK around the phosphorylation site of Thr980:H(p-T)GQ
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	119kD
Uniprot ID	Q9Z2B5
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.