


Product Datasheet

MyoD1/Myf3 Polyclonal Antibody (orb1089049)

Description	MyoD1/Myf3 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse MyoD1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	36 kD
Uniprot ID	P10085
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse & Rat Mus...