

Product Datasheet

Sphingosine 1-Phosphate Lyase 1 antibody (orb108905)

Description

Rabbit polyclonal to Sphingosine 1-Phosphate

Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide derived from the sphingosine 1-phosphate lyase 1 protein.
Preservatives	Provided as ligand affinity purified antibody in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as ligand affinity purified antibody in phosphate buffered saline with 0.08% sodium azide
Concentration	Lot specific, see vial label
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody can be used for Western blotting (1:400 dilution). Other applications not yet tested. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Uniprot ID	O95470
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.

Blocking Peptide

65.5 kDa



Western blot analysis of MDBK cell lysat...